

研究用試薬

Anti-TBE Virus ELISA (IgG) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2661-9601 G	TBE virus	IgG	Ag-coated microplate wells	96 x 01 (96)

Principles of the test: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human antibodies of the IgG class against TBE virus in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with TBE virus antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies (also IgA and IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol								
1. Microplate wells coated with antigens: 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	STRIPS								
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	CAL 1								
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	CAL 2								
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	CAL 3								
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	POS CONTROL								
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	NEG CONTROL								
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	CONJUGATE								
8. Sample buffer ready for use	light blue	1 x 100 ml	SAMPLE BUFFER								
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	WASH BUFFER 10x								
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	SUBSTRATE								
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	STOP SOLUTION								
12. Test instruction.	---	1 booklet									
13. Quality control certificate	---	1 protocol									
<table border="1" style="width: 100%;"> <tr> <td style="width: 10%;">LOT</td> <td style="width: 40%;">Lot</td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> <td style="width: 30%;">Storage temperature</td> </tr> <tr> <td>IVD</td> <td>In vitro determination</td> <td style="text-align: center;"><input type="checkbox"/></td> <td>Unopened usable until</td> </tr> </table>	LOT	Lot	<input type="checkbox"/>	Storage temperature	IVD	In vitro determination	<input type="checkbox"/>	Unopened usable until			
LOT	Lot	<input type="checkbox"/>	Storage temperature								
IVD	In vitro determination	<input type="checkbox"/>	Unopened usable until								

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C, do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.



Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallization occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip: 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Warning: The controls used have been tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays or indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid skin contact.

Preparation and stability of the serum or plasma samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** in sample buffer. For example: dilute 10 µl serum in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: Calibrators and controls are prediluted and ready for use, do not dilute them.



Incubation

For **semiquantative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation: (1st step) Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Wash reagent wells 3 times with 450 µl working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short reaction times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step) Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Empty the wells. Wash as described above.

Substrate incubation: (3rd step) Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C) (protect from direct sunlight).

Stopping the reaction: Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement: **Photometric measurement** of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using the analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on inquiry.

Automated test performance using other fully automated, open system analysis devices is possible, however, the combination should be validated by the user.



Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1-4 is an example for the **semiquantitative analysis** of 24 patient sample (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example for the **quantitative analysis** of 24 patient sample (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The reagent wells are break off format. Therefore, the number of tests performed can be matched to the number of samples, minimizing reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <0.8:	negative
Ratio ≥0.8 to <1.1:	borderline
Ratio ≥1.1:	positive

In cases of borderline test results, an additional patient sample should be taken 7 days later and re-tested in parallel with the first patient sample. The results of both samples allow proper evaluation of titer changes.

Quantitative: The standard curve from which the concentration of antibodies in the serum samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use „point-to-point“ plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.



If the extinction of a serum sample lies above the value of calibrator 1 (200 RU/ml), the result should be given as ">200 RU/ml". It is recommended that the sample be re-tested at a dilution of 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range of non-infected persons (**cut-off value**) recommended by EUROIMMUN is **20 relative units (RU)/ml**. EUROIMMUN recommends interpreting results as follows:

<16 RU/ml:	negative
≥16 to <22 RU/ml:	borderline
≥22 RU/ml:	positive

Evaluation information: For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another the sample should be retested.

For the interpretation of borderline results an investigation using further tests (e.g. avidity determination of antibody class IgG) can be helpful. Diagnosis can be secured by the determination of the titer change in two serum samples taken at an interval of at least 7 days and analysed in parallel.

For diagnosis, the clinical symptoms of the patient should always be taken into account along with the serological results.

Test characteristics

Calibration: As no international reference serum exists for antibodies against TBE virus, the calibration is performed in relative units (RU).

For every group of tests performed, the extinction values of the calibration sera and the relative units and/or ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The activity of the enzyme used is temperature-dependent and the extinction values may vary if a thermostat is not used. The higher the room temperature (+18°C to +25°C) during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibration sera are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.



Antigen: The reagent wells were coated with inactivated TBE virus antigens of the strain "K23" isolated from TBE virus infected PCEC purified chicken embryo cells.

Linearity: The linearity of the Anti-TBE Virus ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95 . The Anti-TBE Virus ELISA (IgG) is linear at least in the tested concentration range (2 - 200 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The detection limit of the Anti-TBE Virus ELISA (IgG) is 0.7 RU/ml.

Cross reactivity: Cross reactivities to other flaviviruses cannot be excluded. They were recognized with anti-West Nile positive and anti-Dengue virus positive samples.

The quality of the antigen used ensures a high specificity of the ELISA. Sera from patients with infections caused by various agents were investigated with the EUROIMMUN Anti-TBE Virus ELISA (IgG).

Antibodies against	n	Anti-TBE Virus ELISA (IgG)
Adenovirus	12	0%
Chlamydia pneumoniae	10	0%
CMV	12	0%
EBV-CA	12	0%
Helicobacter pylori	12	0%
HSV-1	12	0%
Influenza virus A	12	0%
Influenza virus B	12	0%
Measles virus	12	0%
Mumps virus	12	0%
Mycoplasma pneumoniae	12	0%
Parainfluenza virus Pool	12	0%
RSV	12	0%
Rubella virus	12	0%
Toxoplasma gondii	9	0%
VZV	12	0%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for hemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Serum	Mean value (RU/ml)	CV (%)
1	107	3.6
2	144	2.6
3	182	2.0

<i>Inter-assay variation, n = 4 x 6</i>		
Serum	Mean value (RU/ml)	CV (%)
1	113	5.9
2	150	4.3
3	185	3.5



Specificity and sensitivity: 10 clinically characterized patient samples (Interlaboratory test samples from INSTAND, Germany) were examined with this EUROIMMUN ELISA. The test showed a specificity and a sensitivity of 100%, respectively.

n = 10		INSTAND	
		positive	negative
ELISA	positive	7	0
	negative	0	3

Immunization monitoring: 47 patient samples were investigated in a vaccination study with this EUROIMMUN ELISA. Blood was drawn from all patients before and after vaccination. In 44 patients (93.6%) a seroconversion of class IgG and in 40 patients (85.1%) a seroconversion of class IgM occurred. The total results showed that antibodies against TBE virus could be detected in 46 patients (97.9%) after vaccination.

n = 47	IgG positive	IgM positive	IgG and/or IgM positive
before vaccination (neg. patient)	0	0	0
after vaccination	44	40	46

Reference range: The levels of anti-TBE virus antibodies (IgG) were analyzed with this EUROIMMUN ELISA in a panel of 300 healthy blood donors. With a cut-off of 20 RU/ml, 3.7% of the blood donors were anti-TBE virus positive (IgG) which reflects the known percentage of infections in adults.

Clinical significance

Tick-borne encephalitis (TBE) is a systemic infection caused by the human pathogen TBE virus [1]. The coated, single-stranded RNA virus belongs to the flavivirus family [2]. Three subtypes of TBE virus are currently known: type 1 (Western subtype), type 2 (Siberian subtype) and type 3 (Far Eastern subtype) [3, 4].

The Western subtype, which is found in Europe, is transmitted via bites from infected ticks. The most important vectors are types of the genus *Ixodes*, with *Ixodes ricinus* being the most common type in middle Europe [2]. Further types are *Ixodes persulcatus*, which is predominantly found in East Europe and Russia and is a vector for TBE type 2 virus, and more rarely also the genera *Rhipicephalus*, *Dermacentor*, *Haemaphysalis*, *Amblyomma* and from the family of soft ticks the genera *Argas* and *Ornithodoros* [5].

The greatest risk of tick bites occurs in bushy areas, on the edge of woods or in high grass. This is because small mammals and deer, which are the main hosts and primary reservoirs of the blood-sucking pathogen, live in these areas. There are no ticks at heights above 1000 meters [6]. In high-risk areas the percentage of ticks infected with TBE is around one to five percent [6].

With the tick bite the TBE virus is transferred from the tick's salivary gland into the victim's skin, and from here into nearby lymph nodes. From here it moves into further organs, such as connective tissue, skeletal muscle, myocardium, smooth muscle and also monocytes and phagocytes, where the virus multiplies rapidly (primary viraemia) [2, 7]. In a second viraemic phase TBE virus moves into the central nervous system, from where it can enter the brain.

Transmission from virus-infected milk products can occur, although very rarely, via milk from infected goats and sheep, and in exceptional cases from cows [4, 8]. Transmission from person to person does not occur [6].

In a large proportion of patients no verifiable disease symptoms appear [2]. Around 10%-30% of infected persons show symptoms two to twenty days after infection. These include flu-like symptoms such as fever, head and joint aches, and gastrointestinal complaints such as nausea, vomiting, stomach ache and diarrhoea. These symptoms all clear up within a few days [2].

In around 10% of symptomatic patients a second phase with infection of the CNS occurs around one week after defervescence. Symptoms are fever peaks of up to 40°C and signs of brain and meningeal involvement such as headache, vomiting and meningeal irritation.



If the meningoencephalitis progresses, it can impair consciousness to the point of coma and paralysis [2, 9]. 44% of phase 2 patients develop meningitis, 42% meningoencephalitis and 14% myelitis/radiculitis [2, 8, 10]. These symptoms can last several months [9]. Subtype 2 of the TBE virus generally causes more severe symptoms with a higher mortality than subtype 1 [2, 3].

The number of TBE infections has risen continually in recent years in all affected countries [11]. The frequency of TBE infections occurring per 100,000 inhabitants is currently 0.27 in Germany, 3.4 in Sweden, 6.0 in Austria and Switzerland, 24.6 in the Czech Republic and Slovakia, and 50-100 in parts of Russia and West Siberia, [1, 4, 5, 9, 12, 13, 14, 15, 16, 17, 18]. In Europe the greatest affected counties are the Baltic states Lithuania, Estonia and Latvia, where all three subtypes of TBE are in circulation [16, 19]. Because of the often unspecific flu-like symptoms, there are probably a high number of unidentified infections.

The diagnostic method of choice is the demonstration of TBE-virus-specific IgM and IgG in serum or CSF using ELISA [6, 20]. These antibodies can be detected 7-10 days after infection or at the start of the second stage of illness [21, 22]. In cases without clear CNS symptoms (meningitis, encephalomyelitis, myelitis) detection of IgM alone in blood is considered insufficiently specific, and demonstration of a 4-fold titer increase between two serum samples is recommended [5, 6]. Early diagnosis is also possible with the detection of low-avidity specific IgG antibodies [2]. Alternatively, simultaneous detection of TBE-specific IgM and IgG in serum increases the specificity [6]. The detection of specific IgM and IgG antibodies in CSF using ELISA and IIFT indicates CNS involvement [12, 14, 20, 22, 23, 24]. With positive serological CSF results, the CSF also shows pleocytosis and an increase in protein after the second fever peak [5, 20].

ELISA is the method of choice for monitoring the human immune reaction following TBE virus vaccination [25, 26, 27]. TBE vaccination usually leads to a positive antibody titer [26, 28, 29]. Detectable specific IgG antibodies persist with lifelong immunity. In endemic areas of Europe they are detectable in 14%-42% of the population [14]. Special procedures for direct virus detection are the time-consuming cell culture method, nRT-PCR (nested reverse transcriptase polymerase chain reaction) and Westernblot [6, 17, 30].

When TBE virus is detected directly or indirectly and as long as an acute infection is indicated, the infection is notifiable in Germany. The public health authority also commissions various evaluations itself [6].

There is no causal treatment for TBE. Specific antiviral medications do not exist. If the disease has manifested, only symptomatic treatment to reduce individual symptoms is possible [2]. In severe cases intensive medical treatment is necessary.

Prognosis is generally good, especially in children and young adults [10, 31]. The majority of cases of the illness heal without any complications. However, in 10%-30% of symptomatic cases neurological damage of various degrees remains [8, 32]. This includes paralysis, balance difficulties, epilepsy, hearing difficulties and memory and concentration problems [32, 33]. An infection confers lifelong immunity, including against other subtypes of TBE virus. One to two percent of patients with meningoencephalitis die [4, 8, 32].

Very important in differential diagnostics are encephalitides with symptoms similar to those of TBE infections caused by viruses (e.g. adenovirus, Coxsackie virus, cytomegaly virus, Dengue virus, echovirus, Herpes simplex virus, HIV, Japanese encephalitis virus, mumps virus, Powassan virus, poliomyelitis virus, varicella zoster virus, West Nile virus) and other pathogens (e.g. haemophilus influenzae, neisseria meningitidis, streptococcus pneumoniae, borrelia burgdorferi, leptospira, toxoplasma gondii, trichinella spiralis) [34]. Various tests are available for differential diagnosis, particularly BIOCHIP-Mosaics™ for infectious serology (e.g. CNS PROFILE), indirect immunofluorescence (infectious serology), bacteria, microplate ELISA and EUROLINE-WB [35, 36].

Active immunization is the most important preventative measure [4, 8, 9, 26, 37, 38, 39]. It is strongly recommended for all people who live in or visit high-risk areas, including children, adolescents and elderly people, because new vaccines without protein-containing stabiliser allow a risk-free application [28, 29, 33, 39, 40, 41]. Passive immunization is possible up to 3 days following the tick bite (postexposure immunoprophylaxis) [40]. A parallel active immunization is also recommended prophylactically [6, 40]. General preventative protection measures such as exposure prophylaxis should also not be neglected [2].



Reference list

1. Haglund M, Vene S, Forsgren M, Gunther G, Johansson B, Niedrig M, Plyusnin A, Lindquist L, Lundkvist A. **Characterisation of human tick-borne encephalitis virus from Sweden.** J Med Virol 71 (2003) 610-621.
2. Falke D. **Virus der Frühsommer-Meningoenzephalitis (FSME).** In: Hahn H, Falke D, Kaufmann SHE, Ulmann U. **Medizinische Mikrobiologie und Infektiologie.** Springer Verlag (2004) 527-529.
3. Pogodina VV, Bochkova NG, Karan' LS, Frolova MP, Trukhina AG, Malenko GV, Levina LS, Platonov AE. **Comparative analysis of virulence of the Siberian and Far-East subtypes of the tick-borne encephalitis virus.** [Article in Russian] Vopr Virusol 49 (2004) 24-30.
4. Heinz FX, Kunz C. **Tick-borne encephalitis and the impact of vaccination.** Arch Virol Suppl 18 (2004) 201-205.
5. Charrel RN, Attoui H, Clegg JC, Deubel V, Frolova TV, Gould EA, Gritsun TS, Heinz FX, Labuda M, Lashkevich VA, Loktev V, Lundkvist A, Lvov DV, Mandl CW, Niedrig M, Papa A, Petrov VS, Plyusnin A, Randolph S, Suss J, Zlobin VI, de Lamballerie X. **Tick-borne virus diseases of human interest in Europe.** Clin Microbiol Infect 10 (2004) 1040-1055.
6. Robert Koch Institut. **Frühsommer-Meningoenzephalitis (FSME).** RKI-Ratgeber Infektionskrankheiten – Merkblätter für Ärzte (2004).
7. Gelpi E, Preusser M, Laggner U, Garzuly F, Holzmann H, Heinz FX, Budka H. **Inflammatory response in human tick-borne encephalitis: analysis of postmortem brain tissue.** J Neurovirol 12 (2006) 322-327.
8. Lademann M, Wild B, Reisinger EC. **Tick-borne encephalitis (FSME)--how great is the danger really?** [Article in German] MMW Fortschr Med 10 (2003) 45, 47-99.
9. Dumpis U, Crook D, Oksi J.: **Tick-borne encephalitis (review).** Clin Inf Dis 28 (1999) 882–890.
10. Kunze U, Asokliene L, Busse A, Chmelik V, Heinz FX, Hingst V, Kadar F, Kaiser R, Kimmig P, Kraigher A, Krech T, Linnquist L, Rosenfeldt V, Ruscio M, Sandell B, Salzer H, Strle F, Suss J, Zilmer K, Mutz I. **Tick-borne encephalitis in childhood--consensus 2004.** Wien Med Wochenschr 154 (2004) 242-245.
11. Kreuzkamp B. **Early-summer meningoencephalitis cases increase in Europe. Vaccination time is now!** [Article in German] MMW Fortschr Med 148 (2006) 17.
12. Kaiser R, Kern A, Kampa D, Neumann-Haefelin D. **Prevalence of antibodies to Borrelia burgdorferi and tick-borne encephalitis virus in an endemic region in southern Germany.** Zentralbl Bakteriologie 286 (1997) 534-541.
13. Baumberger P, Krech T, Frauchiger B. **Development of early-summer meningoencephalitis (FSME) in the Thurgau region 1990-1995--a new endemic area?** [Article in German] Schweiz Med Wochenschr 126 (1996) 2072-2077.
14. Treib J, Woessner R, Grauer MT, Mueller-Reiland D, Haass A, Schimrigk K. **Prevalence of antibodies to tick-borne encephalitis virus and Borrelia burgdorferi sensu lato in samples from patients with abnormalities in the cerebrospinal fluid.** Zentralbl Bakteriologie 288 (1998) 253-266.
15. Labuda M, Eleckova E, Lickova M, Sabo A. **Tick-borne encephalitis virus foci in Slovakia.** Int J Med Microbiol 291 (2002) 43-47.
16. Suss J, Schrader C, Abel U, Bormane A, Duks A, Kalnina V. **Characterization of tick-borne encephalitis (TBE) foci in Germany and Latvia (1997-2000).** Int J Med Microbiol 291 (2002) 34-42.
17. Suss J, Klaus C, Diller R, Schrader C, Wohanka N, Abel U. **TBE incidence versus virus prevalence and increased prevalence of the TBE virus in Ixodes ricinus removed from humans.** Int J Med Microbiol 296 (2006) 63-68.



18. Rizzoli A, Rosa R, Mantelli B, Pecchioli E, Hauffe H, Tagliapietra V, Beninati T, Neteler M, Genchi C. **Ixodes ricinus, transmitted diseases and reservoirs.** [Article in Italian] *Parassitologia* 46 (2004) 119-122.
19. Bormane A, Lucenko I, Mavtchoutko V, Ranka R, Salmina K, Baumanis V. **Vectors of tick-borne diseases and epidemiological situation in Latvia in 1993-2002.** *Int J Med Microbiol* 293 (2004) 36-47.
20. Kaiser R, Holzmann H. **Laboratory findings in tick-borne encephalitis--correlation with clinical outcome.** *Infection* 28 (2000) 78-84.
21. Niedrig M, Avsic T, Aberle SW, Ferenczi E, Labuda M, Rozentale B, Donoso Mantke O. **Quality control assessment for the serological diagnosis of tick borne encephalitis virus infections.** *J Clin Virol* 2007 Jan 29 [Epub ahead of print]
22. Gunther G, Haglund M, Lindquist L, Skoldenberg B, Forsgren M. **Intrathecal IgM, IgA and IgG antibody response in tick-borne encephalitis. Long-term follow-up related to clinical course and outcome.** *Clin Diagn Virol* 8 (1997) 17-29.
23. Sonnenberg* K, Niedrig M, Steinhagen* K, Rohwader* E, Meyer* W, Schlumberger* W, Müller-Kunert* E, Stöcker* W (*EUROIMMUN AG). **State-of-the-art serological techniques for detection of antibodies against tick-borne encephalitis virus.** *Int J Med Microbiol* 293 (2004) 148-151.
24. EUROIMMUN AG. Stöcker W, Müller M. **Lichtquelle für ein Auflichtfluoreszenzmikroskop.** EUROIMMUN Medizinische Labordiagnostika AG. Eingetragenes deutsches Gebrauchsmuster DE202004010121U (2004).
25. Vene S, Haglund M, Lundkvist A, Lindquist L, Forsgren M. **Study of the serological response after vaccination against tick-borne encephalitis in Sweden.** *Vaccine* 25 (2007) 366-372.
26. Rendi-Wagner P, Kundi M, Zent O, Dvorak G, Jaehrig P, Holzmann H, Mikolasek A, Kollaritsch H. **Persistence of protective immunity following vaccination against tick-borne encephalitis--longer than expected?** *Vaccine* 22 (2004) 2743-2749.
27. Loew-Baselli A, Konior R, Pavlova BG, Fritsch S, Poellabauer E, Maritsch F, Harmacek P, Krammer M, Barrett PN, Ehrlich HJ; FSME-IMMUN study group. **Safety and immunogenicity of the modified adult tick-borne encephalitis vaccine FSME-IMMUN: results of two large phase 3 clinical studies.** *Vaccine* 24 (2006) 5256-5263.
28. Ehrlich HJ, Pavlova BG, Fritsch S, Poellabauer EM, Loew-Baselli A, Obermann-Slupetzky O, Maritsch F, Cil I, Dorner F, Barrett PN. **Randomized, phase II dose-finding studies of a modified tick-borne encephalitis vaccine: evaluation of safety and immunogenicity.** *Vaccine* 22 (2003) 217-223.
29. Kunz C. **TBE vaccination and the Austrian experience.** *Vaccine* 21 (2003) 50-55.
30. Suss J, Schrader C, Falk U, Wohanka N. **Tick-borne encephalitis (TBE) in Germany--epidemiological data, development of risk areas and virus prevalence in field-collected ticks and in ticks removed from humans.** *Int J Med Microbiol* 293 (2004) 69-79.
31. Kaiser R. **Frühsommer-Meningoenzephalitis. Prognose für Kinder und Jugendliche günstiger als für Erwachsene.** *Deutsches Ärzteblatt* 101 (2004) 1822-1826.
32. Gelpi E, Preusser M, Garzuly F, Holzmann H, Heinz FX, Budka H. **Visualization of Central European tick-borne encephalitis infection in fatal human cases.** *J Neuropathol Exp Neurol* 64 (2005) 506-512.
33. Gleixner, Müller, Wirth. **Neurologie und Psychiatrie.** ISBN 3929851539, 4. Aufl. (2004/5) 116.
34. Rendi-Wagner P. **Risk and Prevention of Tick-borne Encephalitis.** *J Travel Med* 11 (2004) 307-312.
35. EUROIMMUN AG. Stöcker W, Teegen B, Meyer W, Müller-Kunert E, Proost S, Schlumberger W, Sonnenberg K. **Differenzierte Autoantikörper-Diagnostik mit BIOCHIP-Mosaiken.** In: Conrad, K. (Hrsg.): **Autoantikörper.** Pabst-Verlag (1998) 78-99.

36. Sonnenberg* K, Hinrichs E, Müller-Kunert* E, Schlumberger* W, Stöcker* W (*EUROIMMUN AG). **Einsatz der BIOCHIP-Technologie in der Serologie von Influenza-, EBV- und Legionella-Infektionen.** Poster zum 3. Deutschen Kongress für Infektions- und Tropenmedizin, 15. bis 18. März 1995 in Berlin-Hohenschönhausen. Abstrakt im Kongressband (1995).
37. Barrett PN., Dorner F, Plotkin SA. **Tick-borne encephalitis vaccine.** In: Plotkin SA, Orenstein (Eds). Vaccines. W.B. Saunders Company (1999).
38. Rendi-Wagner P, Kundi M, Zent O, Banzhoff A, Jaehnig P, Stemberger R, Dvorak G, Grumbeck E, Laaber B, Kollaritsch H. **Immunogenicity and safety of a booster vaccination against tick-borne encephalitis more than 3 years following the last immunisation.** Vaccine 23 (2004) 427-434.
39. Zenz W, Pansi H, Zoehrer B, Mutz I, Holzmann H, Kraigher A, Berghold A, Spork D. **Tick-borne encephalitis in children in Styria and Slovenia between 1980 and 2003.** Pediatr Infect Dis J 24 (2005) 892-896.
40. Robert Koch Institut. Empfehlungen der Ständigen Impfkommission (STIKO) am Robert Koch-Institut. **FSME.** Robert Koch-Institut. Epidemiologisches Bulletin 30 (2005).
41. Zent O, Hennig R, Banzhoff A, Bröker M. **Protection against Tick-Borne Encephalitis with a New Vaccine Formulation Free of Protein-Derived Stabilizers.** Journal of Travel Medicine 12 (2005) 85-93.