



Short instructions for the indirect immunofluorescence test

FA 112I-####-51

Anti-GABA A Receptor (GABARA1/B3) IIFT

ORDER NO.	ANTIBODIES AGAINST	SUBSTRATE	SPECIES	FORMAT SLIDES x FIELDS
FA 112I-1005-51 FA 112I-1010-51	GABA A receptor (GABARA1/B3) ---	Transfected cells Control transfection	EU 90	10 x 05 (050) 10 x 10 (100)



Fig.: Schematic illustration of a slide to demonstrate the arrangement of the BIOCHIPs

Recommended sample dilution

BIOCHIPs	Substrate	Recommended dilution
		IgG
1	GABA A receptor (GABARA1/B3)	1:10 (serum, plasma) or undiluted (cerebrospinal fluid)
2	EU 90 (control transfection)	

Fluorescence pattern (positive reaction):

Autoantibodies against **GABA A receptor (GABARA1/B3)** react specifically with the corresponding transfected cells. They induce a smooth to fine-speckled fluorescence of the cell with an accentuation of the cell membrane, while the area of the cell nucleus is only slightly stained.

Some cells are not transfected and show no specific fluorescence.

If all cells are stained, i.e. also control-transfected cells, antibodies against other cell antigens are present.

For further information refer to the instructions for use FA_0002_A_UK_BXX.

For research use only