

研究用

IIFT: Autoimmune Encephalitis Mosaic 6

ORDER NO.	ANTIBODIES AGAINST	SUBSTRATE	SPECIES	
FA 112d-6	glutamate receptor (type NMDA) glutamate receptor (type AMPA1/2) contactin-associated protein 2 (CASPR2) dipeptidyl aminopeptidase-like protein 6 (DPPX) leucine-rich glioma-inactivated protein 1 (LGI1) GABA _B receptor (GABAR1/B2)	transfected cells	EU 90	10 x 05 (050) 10 x 10 (100)

Indication: This test kit provides qualitative or semiquantitative in vitro determination of human antibodies of immunoglobulin class IgG against neuronal antigens in patient samples to support the diagnosis of neurological diseases (encephalitis). The fluorescence is either evaluated using the fluorescence microscope (specifications see chapter "Incubation", section "Evaluation") or, following automated image recording by the EUROPattern microscope at the computer screen, optionally supported by the EUROPattern Classifier software. The product is designed for use as **VD**. If cerebrospinal fluid is used as the sample material, the results for the parameters dipeptidyl aminopeptidase-like protein 6 (DPPX), leucine-rich glioma-inactivated protein 1 (LGI1) and glutamate receptor (type AMPA1/2) must be confirmed using a CE-labelled test.

Method comparison – Specificity and sensitivity:

Anti-Glutamate receptor (type NMDA) IIFT (*serum/plasma*)

Overview of the tested samples / reference test system:	n
Samples of patients with anti-NMDA receptor encephalitis (n = 29) and a control group (16 patients with other encephalopathies, 1 patient with cerebellar degeneration, 1 patient with retinopathy; origin: Prof. J. Dalmau, USA). The samples were serologically precharacterised as positive or negative by in-house cell-based assay*.	47
Number of samples	47

* Literature reference: Prüss H, Dalmau J, Harms L, Holtje M, Ahnert-Hilger G, Borowski K, Stoecker W, Wandinger KP. **Retrospective analysis of NMDA receptor antibodies in encephalitis of unknown origin.** *Neurology* 75 (2010) 1735-1739.

n = 47		In-house Anti-NMDA-R CBA	
		positive	negative
EUROIMMUN Anti-Glutamate receptor (type NMDA) IIFT	positive	29	0
	negative	0	18

Specificity	100%
Sensitivity	100%

Anti-GABAR β 1/ β 2 IIFT (serum/plasma)

Overview of the tested samples / reference test system:	n
Samples of patients with limbic encephalitis, positive serological precharacterisation by in-house methods (origin: Prof. Dalmau, Hospital of the University of Philadelphia, Dept. of Neurology, Philadelphia, USA).	17
Number of samples	17

n = 17		Reference	
		positive	negative
EUROIMMUN Anti-GABARβ1/β2 IIFT	positive	14	0
	negative	3	0

Specificity	---
Sensitivity	82.4%

Anti-GABAR β 1/ β 2 IIFT (CSF)

Overview of the tested samples / reference test system:	n
Samples of patients with limbic encephalitis, antibody positive precharacterisation by in-house methods (origin: Prof. Dalmau Hospital of the University of Philadelphia, Dept. of Neurology, Philadelphia, USA). 6/8 of the associated serum samples were anti-GABA β receptor positive, no serum material was available for two CSF samples.	10
Number of samples	10

n = 10		Reference	
		positive	negative
EUROIMMUN Anti-GABARβ1/β2 IIFT	positive	6	0
	negative	4	0

Specificity	---
Sensitivity	60%

Anti-LGI1 IIFT (serum/plasma)

Overview of the tested samples / reference test system:	n
Samples of patients with positive or negative serological precharacterisation. The precharacterisation was carried out by means of different in-house methods (in-house-CBA and tissue section) (origin: Prof. Angela Vincent, University of Oxford, United Kingdom).	21
Number of samples	21

n = 21		In-house IFT	
		positive	negative
EUROIMMUN Anti-LGI1 IIFT	positive	8	0
	negative	0	13

Specificity	100%
Sensitivity	100%

Anti-CASPR2 IIFT (serum/plasma)

Overview of the tested samples / reference test system:		n
1. Preselected samples of patients with limbic encephalitis, which were serologically precharacterised by in-house IFT (origin: Prof. Angela Vincent, University of Oxford, United Kingdom).		9
2. Samples of patients with various neurological disorders (incl. limbic encephalitis, epilepsy) and healthy controls, which were serologically precharacterised by in-house IFT (origin: Prof. Angela Vincent, University of Oxford, United Kingdom).		22
Number of samples		31

n = 31		In-house IFT	
		positive	negative
EUROIMMUN Anti-CASPR2 IIFT	positive	14	0
	negative	0	17

Specificity	100%
Sensitivity	100%

Anti-Glutamate receptor (type AMPA1/2) IIFT (serum/plasma)

Overview of the tested samples / reference test system:		n
Preselected samples of patients with positive or negative serological precharacterisation. The precharacterisation was carried out using an in-house IFT (origin: Prof. Josep Dalmau, Hospital Clinic, Dept. of Neurology, Center of Neuroimmunology and Paraneoplastic Disorders, Barcelona, Spain).		13
Number of samples		13

n = 13		In-house IFT	
		positive	negative
EUROIMMUN Anti-AMPA1/R2 IIFT	positive	5	0
	negative	0	8

Specificity	100%
Sensitivity	100%

Anti-DPPX IIFT (serum/plasma)

Overview of the tested samples / reference test system:		n
1. Samples of patients with positive or negative serological precharacterisation were investigated. The precharacterisation was done using an in-house method (origin: Prof. Josep Dalmau, Hospital Clinic, Dept. of Neurology, Center of Neuroimmunology and Paraneoplastic Disorders, Barcelona Spain).		8
2. Samples of patients with positive serological precharacterisation were investigated. The precharacterisation was done using an IFT (tissue sections; origin: Prof. Vanda Lennon, Department of Laboratory Medicine and Pathology, Neuroimmunology Laboratory, Mayo Clinic, USA).		3
Number of samples		11

n = 11		Reference	
		positive	negative
EUROIMMUN Anti-DPPX IIFT	positive	8	0
	negative	0	3

Specificity	100%
Sensitivity	100%

Literature

- Bettler B, Kaupmann K, Mosbacher J, Gassmann M. **Molecular structure and physiological functions of GABA(B) receptors**. *Physiol Rev* 84 (2004) 835-867.
- Boronat A, Gelfand JM, Gresa-Arribas N, Jeong HY, Walsh M, Roberts K, Martinez-Hernandez E, Rosenfeld MR, Balice-Gordon R, Graus F, Rudy B, Dalmau J. **Encephalitis and antibodies to dipeptidyl-peptidase-like protein-6, a subunit of Kv4.2 potassium channels**. *Ann Neurol* 73 (2013) 120-128.
- EUROIMMUN AG. Saschenbrecker S, Rentzsch K, Probst C, Komorowski L, Stoecker W. **Autoantikörperdiagnostik in der Neurologie mittels nativer und rekombinanter Antigen-substrate**. *Der Nervenarzt* 84 (2013) 471-476.
- EUROIMMUN AG, Wandinger KP. **Klinik und Labordiagnostik paraneoplastischer neurologischer Syndrome**. *MTA Dialog* 4 (2008) 266-270.
- EUROIMMUN AG. Stöcker W, Schlumberger W. **Alle Beiträge zum Thema Autoimmun-diagnostik**. In: Gressner A, Arndt T. (Hrsg.). *Lexikon der Medizinischen Laboratoriumsdiagnostik*. 2. Aufl., Springer Medizin Verlag, Heidelberg (2012).
- Finke C, Kopp UA, Prüss H, Dalmau J, Wandinger* KP, Ploner CJ. (*EUROIMMUN AG). **Cognitive deficits following anti-NMDA receptor encephalitis**. *J Neurol Neurosurg Psychiatry* 83 (2012) 195-198.
- Finke C, Prüss H, Scheel M, Ostendorf F, Harms L, Borowski* K, Wandinger* KP, Ploner CJ. (*EUROIMMUN AG). **Anti-NMDA receptor antibodies in a case of MELAS syndrome**. *J Neurol* 259 (2012) 582-584.
- Fukata Y, Lovero KL, Iwanaga T, Watanabe A, Yokoi N, Tabuchi K, Shigemoto R, Nicoll RA, Fukata M. **Disruption of LGI1-linked synaptic complex causes abnormal synaptic transmission and epilepsy**. *Proc Natl Acad Sci U.S.A.* 107 (2010) 3799-3804.
- Kim J, Nadal MS, Clemens AM, Baron M, Jung SC, Misumi Y, Rudy B, Hoffman DA. **Kv4 accessory protein DPPX (DPP6) is a critical regulator of membrane excitability in hippocampal CA1 pyramidal neurons**. *J Neurophysiol* 100 (2008) 1835-1847.
- Kumar MA, Jain A, Dechant VE, Saito T, Rafael T, Aizawa H, Dysart KC, Katayama T, Ito Y, Araki N, Abe T, Balice-Gordon R, Dalmau J. **Anti-N-methyl-D-aspartate receptor encephalitis during pregnancy**. *Arch Neurol* 67 (2010) 884-887.

- Lai M, Hughes EG, Peng X, Zhou L, Gleichman AJ, Shu H, Mata S, Kremens D, Vitaliani R, Geschwind MD, Bataller L, Kalb RG, Davis R, Graus F, Lynch DR, Balice-Gordon R, Dalmau J. **AMPA receptor antibodies in limbic encephalitis alter synaptic receptor location.** Ann Neurol 65 (2009) 424-434.
- Lancaster E, Lai M, Peng X, Hughes E, Constantinescu R, Raizer J, Friedman D, Skeen MB, Grisold W, Kimura A, Ohta K, Iizuka T, Guzman M, Graus F, Moss SJ, Balice-Gordon R, Dalmau J. **Antibodies to the GABA(B) receptor in limbic encephalitis with seizures: case series and characterisation of the antigen.** Lancet Neurol 9 (2010) 67-76.
- Misawa T, Mizusawa H. **Anti-VGKC antibody-associated limbic encephalitis/Morvan syndrome.** Brain Nerve 62 (2010) 339-345.
- Mombereau C, Kaupmann K, Froestl W, Sansig G, Van der PH, Cryan JF. **Genetic and pharmacological evidence of a role for GABA(B) receptors in the modulation of anxiety- and antidepressant-like behavior.** Neuropsychopharmacology 29 (2004) 1050-1062.
- Nobile C, Michelucci R, Andreazza S, Pasini E, Tosatto SC, Striano P. **LGI1 mutations in autosomal dominant and sporadic lateral temporal epilepsy.** Hum Mutat 30 (2009) 530-536.
- Prüss H, Dalmau J, Harms L, Höltje M, Ahnert-Hilger G, Borowski* K, Stoecker* W, Wandinger* KP. (*EUROIMMUN AG). **Retrospective analysis of NMDA receptor antibodies in encephalitis of unknown origin.** Neurology 75 (2010) 1735-1739.
- Prüss H, Höltje M, Maier N, Gomez A, Buchert R, Harms L, Ahnert-Hilger G, Schmitz D, Terborg C, Kopp U, Klingbeil* C, Probst* C, Kohler S, Schwab JM, Stöcker* W, Dalmau J, Wandinger* KP. (*EUROIMMUN AG). **IgA NMDA receptor antibodies are markers of synaptic immunity in slow cognitive impairment.** Neurology 78 (2012) 1743-1753.
- Seikel E, Trimmer JS. **Convergent modulation of Kv4.2 channel alpha subunits by structurally distinct DPPX and KChIP auxiliary subunits.** Biochemistry 48 (2009) 5721-5730.
- Sirerol-Piquer MS, Ayerdi-Izquierdo A, Morante-Redolat JM, Herranz-Perez V, Favell K, Barker PA, Perez-Tur J. **The epilepsy gene LGI1 encodes a secreted glycoprotein that binds to the cell surface.** Hum Mol Genet 15 (2006) 3436-3445.
- Tan KM, Lennon VA, Klein CJ, Boeve BF, Pittock SJ. **Clinical spectrum of voltage-gated potassium channel autoimmunity.** Neurology 70 (2008) 1883-1890.
- Tuzun E, Dalmau J. **Limbic encephalitis and variants: classification, diagnosis and treatment.** Neurologist. 13 (2007) 261-271.
- Vincent A, Irani SR. **Caspr2 antibodies in patients with thymomas.** J Thorac Oncol 5 (2010) S277-S280.
- Wandinger* KP, Dalmau J, Borowski* K, Probst* C, Fechner* K, Stöcker* W. (*EUROIMMUN AG). **Recombinant immunofluorescence assay for the detection of anti-glutamate receptor (type NMDA) antibodies in the differential diagnosis of autoimmune encephalopathies.** In: Conrad K et al. (Hrsg.). From Pathogenesis to Therapy of Autoimmune diseases: Autoantigens, Autoantibodies, Autoimmunity. Pabst Science Publishers (2009) 434-435.
- Wandinger* KP, Klingbeil C*, Gneiss C, Waters P, Dalmau J, Saschenbrecker* S, Borowski K*, Deisenhammer F, Vincent A, Probst* C, Stoecker* W. (*EUROIMMUN AG). **Neue serologische Marker zur Differentialdiagnose der Autoimmun-Enzephalitis/New serological markers for the differential diagnosis of autoimmune limbic encephalitis.** J Lab Med 35 (2011) 329-342.
- Wandinger* KP, Saschenbrecker* S, Stoecker* W, Dalmau J. (*EUROIMMUN AG). **Anti-NMDA-receptor encephalitis: A severe, multistage, treatable disorder presenting with psychosis.** J Neuroimmunol 231 (2011) 86-91.