




# Immunoblot PreQ

Pre-equipped (PreQ) incubation trays for fully automated processing of EUROLINE blot strips



Autoimmune diagnostics




Infection diagnostics

- Lowest manual workload owing to pre-equipment of incubation trays with EUROLINE blot strips
- Reliable fixation of EUROLINE blot strips
- User-friendly, patented\* click-in system allows flexible combination of the PreQ incubation trays
- Fully automated incubation and evaluation by EUROBlotOne and EUROLineScan software



The **Immunoblot PreQ products** offered by EUROIMMUN are a time-saving and convenient alternative to the manual positioning of EUROLINE blot strips in the incubation tray for use with the fully automated system EUROBlotOne. The individual PreQ incubation trays are pre-equipped with EUROLINE blot strips and can be flexibly combined via a simple click-in system (at least 4 per run). The incubation, drying and subsequent image recording are carried out fully automatically by the EUROBlotOne. The subsequent evaluation and archiving of images is performed reliably and conveniently by the EUROLineScan software.

The optimal drying time for each format is transmitted to the EUROBlotOne and adjusted fully automatically by selecting "Use Immunoblot PreQ" in the EUROLineScan software.



## Currently available Immunoblot PreQ products for autoimmune and infection diagnostics

Order no.	Product
DL 1111-5001-7 G	EUROLINE Paraneoplastic Neurological Syndromes 12 Ag (IgG)
DL 1300-5001-4 G	EUROLINE Autoimmune Liver Diseases (IgG)
DL 1300-5001-9 G	EUROLINE Autoimmune Liver Diseases 9 Ag plus F-Actin (IgG)
DL 1360-0510 A	
DL 1360-5001 G	EUROLINE Autoimmune Gastrointestinal Diseases (IgA, IgG)
DL 1530-5001-4 G	EUROLINE Autoimmune Inflammatory Myopathies 16 (IgG)
DL 1530-5001-7 G	EUROLINE Autoimmune Inflammatory Myopathies 16 Ag et cN-1A (IgG)
DL 1530-5001-9 G	EUROLINE Autoimmune Inflammatory Myopathies 20 Ag(IgG)
DL 1532-5001 G	EUROLINE Systemic Sclerosis (Nucleoli) Profile (IgG)
DL 1532-5001-2 G	EUROLINE Systemic Sclerosis (Nucleoli) Profile 2 (IgG)
DL 1590-5001-3 G	EUROLINE ANA Profile 3 (IgG)
DL 1590-5001-8 G	EUROLINE ANA Profile 1 (IgG)
DL 1590-5001-23 G	EUROLINE ANA Profile 23 (IgG)
DL 1590-5001-30 G	EUROLINE ANA Profile 3 plus DFS70 (IgG)
DL 1590-5001-32 G	EUROLINE Anti-dsDNA, -Nucleosomes, -Histones, -DFS70 (IgG)
DL 1590-5001-33 G	EUROLINE ANA Profile et Mi-2, Ku, DFS70 (IgG)
DL 1590-5001-35 G	EUROLINE Cytoplasm Profile (IgG)
DL 1590-5001-37 G	EUROLINE Cytoplasm Profile 13 Ag (IgG)
DL 159z-5001 G	EUROLINE Anti-DFS70 (IgG)
DN 2131-5001 G, M	Anti-Borrelia EUROLINE-RN-AT (IgG, IgM)
DN 2131-5001-2 M	Anti-Borrelia EUROLINE-RN-AT-adv (IgM)
DY 2131-4801-1 G, M	Anti-Borrelia EUROLINE-WB (IgG, IgM)
DY 2190-0148-1 A	
DY 2190-4801-1 G	Anti-Chlamydia HP EUROLINE-WB (IgA, IgG)
DY 2321-4801-2 G	EUROLINE-WB Anti-Echinococcus (IgG)
DN 2410-5001-4 G, M	EUROLINE Anti-TO.R.C.H. Profile (IgG, IgM)
DN 2410-5001-11 G	EUROLINE Anti-TO.R.C.H. Profile 10 (IgG)
DN 2525-0510 A	
DN 2525-5001 G	EUROLINE Anti-Hepatitis E Virus (IgA, IgG)
DN 2570-5001 G, M	EUROLINE Anti-Cytomegalovirus (IgG, IgM)
DY 2590-4801 G	Anti-Rubella Virus WESTERNBLOT (IgG)
DN 2606-5001-1 G	EUROLINE Anti-SARS-CoV-2 Profile (IgG)
DN 2790-0510-2 M	
DN 2790-5001-2 G	EUROLINE Anti-EBV Profile 2 (IgM, IgG)
ZG 015o-0101-146	Positioning aid for EUROBlotOne (one-time setup is a prerequisite for using the Immunoblot PreQ)