

euroimmun 

From Revvity

Your partner
for modern allergy
diagnostics.




Increased efficiency

High throughput



- True walk-away system
- 120 sample positions
- 80 results/ hour


True random access



- Random access for all tests of the different diagnostic areas
- Continuous loading

Multi-Discipline Automated System

Efficient and reliable laboratory diagnostics



Easy-to-use instrument

Reduced set-up time



IDS

5 minutes

Common system

> 20 minutes

- Automatic start-up and shut-down
- First test result after 40 minutes

Minimal hands-on time



IDS

15 minutes

Common system

Several hours

- Automatic self-maintenance
- Rapid user maintenance: daily, weekly and monthly
- Ready-to-use reagents

Reagent module

- On-board refrigeration of ready to use immunoassays
- Continuous loading and automated barcode recognition

Efficiency

- Fully automated walk-away system
- Compact, benchtop design
- Efficient random access loading

ids isys

Allergy diagnostics

Efficient and complete solutions for your testing needs



Weed pollens



Fish & shellfish



Microorganisms



Eggs



Fruits



Milk



Epidermals



Grass pollens



Mites



Tree pollens



Legumes, nuts & seeds



Occupational



Insect venoms

- Automated quantitative determination of total and specific IgE, the most important serological markers
- Combine assays from different diagnostic areas (allergy, autoimmunity, therapeutic drug monitoring, endocrinology, alzheimer's disease, infection) in one run to fully utilise existing capacities
- Detection of antibodies against the most common and relevant allergens
- Specific and sensitive tests based on allergen components and extracts
- Development based on recognised guidelines¹
- High correlation with competitor systems

¹. Matricardi PM, Kleine-Tebbe J, Hoffmann HJ, Valenta R, Hilger C, et al. EAACI Molecular Allergology User's Guide. *Pediatr Allergy Immunol.* 2016 May;27 Suppl 23:1-250

Find out more at www.euroimmun.com
or contact us directly:
www.euroimmun.com/contact



euroimmun 

From Revvity

EUROIMMUN Medizinische Labordiagnostika AG
Seekamp 31 · 23560 Lübeck (Germany) · Phone: +49 451 2032-0
www.euroimmun.com

Regulatory status of the products must be verified for the user's individual jurisdiction.
Please contact your country representative for product availability and information.