




# EUROArray STI



## Detection of up to 11 relevant sexually transmitted pathogens in one reaction




BACTERIAL PATHOGENS	PROTOZOON	VIRAL PATHOGENS
<i>Chlamydia trachomatis</i> (CT)	<i>Trichomonas vaginalis</i> (TV)	Herpes simplex virus 1 (HSV-1)
<i>Haemophilus ducreyi</i> (HD)		Herpes simplex virus 2 (HSV-2)
<i>Mycoplasma genitalium</i> (MG)		
<i>Mycoplasma hominis</i> (MH)		
<i>Neisseria gonorrhoeae</i> (NG)		
<i>Treponema pallidum</i> (TP)		
<i>Ureaplasma parvum</i> (UP)		
<i>Ureaplasma urealyticum</i> (UU)		

- Reliable DNA-based detection of 11 sexually transmitted pathogens
- Identification of pathogens which are difficult or impossible to cultivate
- Easy detection of multiple infections

**EUROArray STI**  
(order no. MN 2830-####)



####: For detailed information about available test kit formats visit [www.euroimmun.com](http://www.euroimmun.com): 

- Fully automated evaluation
- Integrated controls for all work steps
- Simple procedure – just one PCR, ready-to-use reagents
- Modular and efficient automation

*Regulatory status of the products must be verified for the user's individual jurisdiction. Please contact your country representative for product availability and information.*